

MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.

VÖLKER

Völker accessories








for S 964, S 962-2, S 962-2 W and S 962-2 Vis-a-Vis hospital beds





















Völker accessories for S 964, S 962-2, S 962-2 W and S 962-2 Vis-a-Vis hospital beds









All Völker accessories attach through simply inserting them into or suspending them from the integrated sockets, holders or rails.

The following pages provide a detailed overview of which accessories are available for which hospital bed.

Description		Product no.	S 964	S 962-2 S 962-2 W	S 962-2 VaV
	Fold-down head or foot board Only available for a bed width of 90 cm	PZK-1035/4	-	●	-
	Lifting pole with handle Material: powder-coated	PZP-3060	●	●	-
		PZK-936/1	●	●	●
	Lifting pole with handle With one-hand adjustment Material: powder-coated	PZP-3060/2	●	●	-
		PZK-936/3	●	●	●
	Document holder	PZK-1070	●	●	●
	Assist rail spacer 1 piece Material: stainless steel and High-Pressure Laminate (HPL)	PZK-969/2	●	●	●
	Assist rail spacer with perforation 1 piece Material: stainless steel and High-Pressure Laminate (HPL)	PZK-969/3	●	●	●
	Assist rail spacer cover Colour: Grey	PZK-982	●	●	●

Description		Product no.	S 964	S 962-2 S 962-2 W	S 962-2 VaV
	Assist rail spacer holder	PZK-979	•	•	•
	Restraint strap holder 1 pair	Easy to adjust	PZK-972/3	-	•
		Screwed tightly in place	PZK-972/4	-	•
	Restraint strap holder 1 pair	PZK-972/5	•	-	-
	Crutch holder 1 pair	PZK-965	•	•	•
	Crutch holder 1 holder fits 2 crutches	PZK-965/2	•	•	•
	Pendant control holder	PZK-1044/3	-	•	•
	Towel holder	PZK-978/2	•	•	•
	IV holder For lifting pole PZP-3060	PZP-3055/2	•	•	-
	IV holder For lifting pole PZK-936	PZK-937/3	•	•	•

Description		Product no.	S 964	S 962-2 S 962-2 W	S 962-2 VaV
	IV pole With one-hand adjustment Maximum load per hook 2 kg Material: stainless steel	PZK-939_R	•	•	•
	IV pole With one-hand adjustment Maximum load per hook 2 kg Material: stainless steel	PZK-939_S	•	•	•
	IV pole (straight) With one-hand adjustment Maximum load per hook 2 kg Material: stainless steel	PZK-939_U	•	•	•
	Material: stainless steel	Length: 40 cm PZK-1083/6	-	•	•
		Length: 60 cm PZK-1083/8	-	•	•
	ISO rail 1 pair Centre, length: 35 cm Material: stainless steel	PZK-1064 M	•	-	-
	ISO rail 1 pair Foot section, length: 20 cm Material: stainless steel	PZK-1064 F	•	-	-
	ISO rail Attached to the backrest at the head section Material: stainless steel	PZK-1091	-	•	•
	Name plate	PZK-950_R	•	•	•
	Equipotential cable Length: 300 cm	PZK-948	•	•	•

Description		Product no.	S 964	S 962-2 S 962-2 W	S 962-2 VaV
	Rhodon basket Material: chrome	PZK-1041	-	•	•
	Writing surface Material: High-Pressure Laminate (HPL) Colour: Grey	PZK-962/2	•	•	•
	Assist rail pads For bed length: 200 cm	Side rail height: 34 cm Side rail height: 37 cm			
		PZK-981_2.0	•	•	•
		PZK-981EN20			
	General-purpose hook 1 pair	PZK-1043	-	•	•
	General-purpose hook 1 pair For ISO rails PZK-1083/PZK-1064	PZK-943	•	•	•
	Urine bottle basket Mounts to accessory rail Material: chrome	PZK-1040	-	•	•
	Urine bottle basket Mounts in socket Material: chrome	PZK-1040/2	•	•	•
	Urine bottle basket For ISO rail PZK-1083 Material: chrome	PZK-940/3	-	•	•

• = available | - = not available

Additional accessories	Product no.	S 964	S 962-2 S 962-2 W	S 962-2 VaV
Arm rests 1 set	PZK-ARVAV	-	-	●
Foot rests 1 set	PZK-FRVAV	-	-	●
Bed sheets 2 pieces, for all bed sizes	PZK-BSVAV	-	-	●

● = available | - = not available

The availability of certain accessories also depends on the year of production and configuration of your bed.

Would you like to know exactly which accessory versions are suitable for your Völker bed? Please note the serial number of your bed and contact your local representative.



Made in Germany for over 100 years. Dedicated to dignified care.

Guided by the principle, Völker develops and produces healthcare and hospital beds that allow patients to recover as fully as possible and that offer residents exceptional comfort. With their extensive functionality, our beds are versatile enough for any hospital or healthcare environment. Their appearance ensures that every room has a pleasant atmosphere. Our beds combine form and function to meet the needs of patients, residents and caregivers alike.

Original Völker accessories make the perfect addition to our beds. As well as meet-

ing all the relevant safety global requirements.

Our manufacturing procedures in Germany ensures the highest quality and processing standards across all our products. With comprehensive guidance (before and after purchasing), efficient logistics, technical service and advanced training capabilities, we offer our customers an all-inclusive service model. This enhances the care process, reduces follow up costs and increases caregiver and user satisfaction - not to mention giving staff more time to provide dignified care.

Play it safe. With Völker Service.

Völker Service maximises the return on your investment in Völker beds. Our Service team will also take responsibility for a large amount of the work and requirements that arise from the German Medical Devices Act and other legislation.

All of our service technicians are qualified and

experienced testers who are very familiar with the medical devices they inspect. They will test and service your beds meticulously, observing all of the legal requirements.

Any defects are usually fixed immediately and therefore at optimum cost.

Völker Service additionally compiles detailed legal documentation and checks all records to ensure that they are complete.

We would be happy to discuss which Völker Service options would best suit your needs. Simply get in touch.



MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.

Völker reserves the right to make changes to the design, specifications and models without prior notice. Products may vary slightly from the colours shown here. Subject to printing errors. All information is non-binding.

Not all products/options are available in all countries. For further information on our products or services, please contact your local Völker representative or visit our website.

©2018 Völker | All rights reserved

VÖLKER

www.voelker.de



**HEALTHCAREBEDS
AUSTRALIA**

info@healthcarebedsaust.com.au