

MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.

VÖLKER

Bedside cabinets and overbed tables

for hospitals and healthcare environments



Hospital bedside cabinets

Model K 931



Model K 940



Model K 941
Model K 941/2



Model K 941 MK



Model

Dimensions in cm (W x D x H)

Overbed tray in cm (W x D)

Height-adjustable overbed tray in cm*

Weight*

Material

Colour

Equipment

Bed tray

Can be used on both sides

Single bottle holder in door (on both sides)

Triple bottle holder in drawer (on both sides)

Castors, dia. 6.5 cm (2 castors with brake)

Door

Drawer instead of door

Open shelf

Bed pan holder

Mount for telephone, can be positioned in one of two places

Bar top

1 towel holder

Lockable upper drawer

3-sided closed rail, height: 68 mm, 3D



Single bottle holder



Triple bottle holder



	K 931	K 940	K 941	K 941/2	K 941 MK with fridge
	48 x 48 x 73	60 x 48 x 57 (or 73)	60 x 48 x 73	60 x 48 x 86	68 x 56 x 79
	-	55 x 35	55 x 35	55 x 35	55 x 35
	-	64 - 98	64 - 98	64 - 98	64 - 98
	25 kg	35 kg	35 kg	38 kg	48 kg
Aluminium (RAL 9006) and high-pressure laminate (HPL)					
See colour collection					

	-	●	●	●	●
	●	●	●	●	-
	○	○	○	○	-
	○	○	○	○	-
	●	●	●	●	●
	●	●	●	●	-
	○	○	○	○	-
	●	-	●	●	-
	○	○	○	○	○
	○	○	○	○	○
	-	●	-	-	-
	○	○	○	○	○
	○	○	○	○	○
	○	-	○	○	○

● = standard | ○ = optional extra | - = not available | X = optional extra, free of charge | * for standard configuration



Bed pan holder



Mount for telephone



Towel holder



3-sided closed rail

Healthcare bedside cabinets

Model P 920



Model P 930



Model P 935



Model P 941 O.T.



Model P 941



Model

Dimensions in cm (W x D x H)

Overbed tray in cm (W x D)

Height-adjustable overbed tray in cm*

Weight*

Material

Design

Colour

Equipment

Bed tray

Can be used on both sides

Single bottle holder in door (on both sides)

Single bottle holder in the door

Triple bottle holder in the drawer

Castors, dia. 6.5 cm (2 castors with brake)

Door instead of drawer

Drawer

Open shelf

Lockable upper drawer

Overbed table

Model
RT-985



	P 920	P 930	P 935	P 941	P 941 O.T.
	42 x 44 x 61	48 x 48 x 61	48 x 48 x 61	60 x 48 x 73	48 x 48 x 74
	-	-	-	55 x 35	-
	-	-	-	64 - 98	-
	24 kg	24 kg	24 kg	33 kg	28 kg
	Decorative chipboard				
	Panelled design	Round-profile design	Square-corner profile	Round-profile design	Round-profile design
	See colour collection				

-	-	-	●	-
-	-	-	●	●
-	-	-	○	○
○	○	○	-	-
○	○	○	-	-
●	●	●	●	●
X	X	X	●	●
●	●	●	-	-
-	-	-	●	●
○	○	○	○	○

● = standard | ○ = optional extra | - = not available | X = optional extra, free of charge | * for standard configuration



Model	RT-985
Dimensions in cm (W x D x H), without railing	98 x 45 x 74 to 117
Weight	20 kg
Material	White aluminium (RAL 9006) and high-pressure laminate (HPL)
Colour	See colour collection

Equipment

Castor, dia. 5.0 cm*	X
Castor, dia. 6.5 cm*	●
Castor, dia. 7.5 cm*	○
Fully rotating tray	X

* = 2 castors with brakes

Völker bedside cabinets – the perfect addition to any Völker bed



Model P 930



Model K 940

Made in Germany for over 100 years. Dedicated to dignified care.

Guided by the principle, Völker develops and produces healthcare and hospital beds that allow patients to recover as fully as possible and that offer residents exceptional comfort. With their extensive functionality, our beds are versatile enough for any hospital or healthcare environment. Their appearance ensures that every room has a pleasant atmosphere. Our beds combine form and function to meet the needs of patients, residents and caregivers alike.

Original Völker accessories make the perfect addition to our beds. As well as meet-

ing all the relevant safety global requirements.

Our manufacturing procedures in Germany ensures the highest quality and processing standards across all our products. With comprehensive guidance (before and after purchasing), efficient logistics, technical service and advanced training capabilities, we offer our customers an all-inclusive service model. This enhances the care process, reduces follow up costs and increases caregiver and user satisfaction - not to mention giving staff more time to provide dignified care.

Play it safe. With Völker Service.

Völker Service maximises the return on your investment in Völker beds. Our Service team will also take responsibility for a large amount of the work and requirements that arise from the German Medical Devices Act and other legislation.

All of our service technicians are qualified and

experienced testers who are very familiar with the medical devices they inspect. They will test and service your beds meticulously, observing all of the legal requirements.

Any defects are usually fixed immediately and therefore at optimum cost.

Völker Service additionally compiles detailed legal documentation and checks all records to ensure that they are complete.

We would be happy to discuss which Völker Service options would best suit your needs. Simply get in touch.



**MADE IN GERMANY FOR OVER 100 YEARS.
DEDICATED TO DIGNIFIED CARE.**

Völker reserves the right to make changes to the design, specifications and models without prior notice. Products may vary slightly from the colours shown here. Subject to printing errors. All information is non-binding.

Not all products/options are available in all countries. For further information on our products or services, please contact your local Völker representative or visit our website.

©2018 Völker | All rights reserved

VÖLKER

www.voelker.de



**HEALTHCAREBEDS
AUSTRALIA**

info@healthcarebedsaust.com.au